



Diagnostica aus Berlin

diaglobal

Diagnostica

Step by step

Control measuring with ERY QS

Control measuring with ERY QS Site 1




1



Withdraw 10µl control solution with capillary from ERY QS bottle

2



Insert capillary in an ERY 142 cuvette

3



Wash out control solution with micropipetter

4



Screw on turquoise-coloured cap again
Turn cuvette upside down several times
Wait 3 min.

For Compur-Minilab-User:


For the control measuring of ERY QS with the photometer Compur Minilab just withdraw 5 µl blood from ERY QS- bottle and assay as usual according to the working instruction of Compur Minilab device.

5



Switch on photometer
Press Enter
Select test „ERY“ with
arrow key
Confirm with Enter

6



Insert not treated cuvette in
Photometer (blank value)
After the signal tone,
remove cuvette

7



Insert blood
sample (pic. 4)
into photometer

8



Read result and compare
with target value.
Max. tolerance from
target value is 8 %.

For Compur-Minilab-User:

For the control measuring of ERY QS with the photometer Compur Minilab just withdraw 5 µl blood from ERY QS- bottle and assay as usual according to the working instruction of Compur Minilab device.